

An Improved Extraction Strategy for Polar Solvent Reactions

With many reactions today being performed using polar, high boiling point solvents such as DMSO, the post-synthesis work-up is not trivial. With a boiling point of 189 °C, DMSO is among the most difficult solvents to evaporate thus requiring chemists to extract with an immiscible solvent with a lower boiling point. However, with most solvents being soluble in DMSO, water must be added before the addition of a suitable extractant such as DCM. This is the approach many chemists follow for this situation.

The typical extraction device is a glass separatory funnel, which all chemists know how to use. However, there are drawbacks to the separatory funnel.

- » Fragility
- » Emulsion formation
- » Gas/pressure evolution
- » Phase interface recognition
- » Flow rate restrictions
- » Apparatus cleanliness/maintenance

Biotage has addressed these limitations by developing efficient and recyclable polypropylene Phase Separators, Figure 1. These vessels contain a hydrophobic frit which allows dense organic solvents such as DCM to flow while retaining water thus eliminating the glass separatory funnel drawbacks.




Figure 1. ISOLUTE® Phase Separator mounted on a Biotage® Gravity Rack.

The ISOLUTE® Phase Separator eliminates the need to monitor flow rate and the organic/aqueous phase interface. Because the reaction mixture can be mixed with water and DCM in the reaction vessel, there is no need to vent pressure or even be concerned with emulsions.

To show their effectiveness, the phase separators were used to extract the reaction product of a succinic acid and benzylamine reaction in DMSO from several by-products, Figure 2.



Figure 2. Reaction of succinic acid and benzylamine in DMSO using a Biotage® Initiator+ microwave reactor.

Materials and Methods

Synthesis and Work-up

Reactor:	Biotage® Initiator+ with 2–5 mL reactor vial
Solvent:	DMSO
Reagents:	Succinic acid (1 eq), benzylamine (2 eq)
Reaction:	15 minutes at 200 °C
Extraction:	ISOLUTE® Phase separator, 25 mL, with Biotage® Gravity Rack

Evaporation

System:	Biotage® V-10 Touch Rapid Solvent Evaporator
----------------	--

Purification

System:	Biotage® Selekt
Column:	Biotage® Sfär C18 12 gram
Dry Load Vessel:	Biotage® Sfär C18 Samplet®
Load Amount:	~21 mg
Solvent A:	Water, deionized
Solvent B:	Methanol
Flow Rate:	30 mL/min.
Equilibration:	10% B for 2 CV
Gradient:	10% B for 1 CV 10–100% B in 10 CV 100% B for 2 CV
Detection:	UV 200–400 nm

Results and Discussion

To evaluate the DMSO reaction mix complexity, a 0.1 mL aliquot was loaded onto a 12 g Sfär C18 column (prior to extraction). The results indicated two major retained compounds and several by-products, Figure 3.




Figure 3. Reversed-phase flash chromatography of the DMSO crude reaction mix prior to extraction shows two major products (fractions 7 and 12) and many by-products. Fraction 1 is DMSO.

To evaluate the ISOLUTE® phase separator's efficiency, a 0.5 mL aliquot was transferred to a 20 mL scintillation vial. To this vial was added ~5 mL of water and 3 mL of DCM. The contents were mixed and poured into a Phase Separator attached to a Biotage® Gravity Rack.

Within seconds the dense DCM extract was collected into a scintillation vial. The crude reaction vial was rinsed again with both water and DCM, transferred to the Phase Separator, and the DCM portion collected. The entire extraction process required only about 2 minutes.

Since DCM is not a suitable injection solvent for reversed-phase, the extract was evaporated on a Biotage® V-10 Touch followed by re-dissolution in a 1:1 mix of methanol and acetone (1 mL). Purification (0.1 mL of re-dissolved extract) was repeated and the results showed removal of most interfering by-products, Figure 4.




Figure 4. Flash chromatography of the DCM extract (dried and re-dissolved in 1:1 methanol/acetone) using reversed-phase shows removal of many of the reaction by-products.

Conclusion

Reactions conducted in high boiling point solvents pose work-up and evaporation challenges requiring a substantial time investment. Biotage Phase Separators simplify and speed compound extraction with equal, if not better, efficiency compared to glass separatory funnels. After use, phase separators can be dried and recycled with other polypropylene lab waste.

EUROPE

Main Office: +46 18 565900
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 800 446 4752
Order Fax: +1 704 654 4917
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: +82 31 706 8500
Fax: +82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

INDIA

Tel: +91 11 45653772
india@biotage.com

Distributors in other regions
are listed on www.biotage.com

Literature Number: AN166

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.