# Extraction of Barbiturates From Human Urine Using ISOLUTE<sup>®</sup> SLE+ Columns with GC-MS Analysis

## Introduction

This application note describes the extraction of a range of barbiturates from human urine using ISOLUTE SLE+ supported liquid extraction columns followed by GC-MS analysis.



This method describes the use of ISOLUTE SLE+ supported liquid extraction 1 mL sample volume columns to extract a range of barbiturates from human urine. The analysis of these analytes was carried out by GC-MS. This simplified and efficient extraction method has significant analyte recoveries ranging from 103-108% with LOQs of 10 ng/mL and RSDs

Figure 1. Structure of Barbituric acid, the basic structure of all barbiturates

ISOLUTE SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquidliquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation time.

<10%.

## Analytes

Butabarbital, butalbarbital, amobarbital, pentobarbital, secobarbital, hexobarbital, phenobarbital.

## **ISOLUTE SLE+ procedure**

ISOLUTE SLE + 1 mL sample volume column, part number 820-0140-C

Sample Pre-treatment:	To 500 $\mu L$ of urine add 100mM ammonium acetate pH 5 (500 $\mu L,~1{:}1,~v/v).$
Sample Load:	Load pre-treated sample $(1 \text{ mL})$ to column followed by a pulse of vacuum to initiate flow and allow to absorb for five minutes.
Analyte Extraction:	Elute with dichloromethane (2.5 mL). Leave to flow under gravity for 5 minutes, then follow with a further aliquot of dichloromethane (2.5 mL) and allow to flow under gravity for a further five minutes, to complete extraction apply a short pulse of vacuum.
Post extraction:	Evaporate to dryness at room temperature (80 L/min) and reconstitute in ethyl acetate (200 $\mu L).$

# Reagents

Ethyl acetate from Fisher Scientific, Loughborough.

Ammonium acetate from Sigma-Aldrich, Gillingham.

Trimethylphenylammonium hydroxide (TMAH) from Sigma-Aldrich, Gillingham.



# **GC Conditions**

Carrier:	Helium 2 mL min <sup>-1</sup> (constant flow)
Inlet:	Splitless, 150 °C
Injection:	In-port flash alkylation: 1 $\mu L$ sample + 1 $\mu L$ 0.2M TMAH (trimethylphenylammonium hydroxide) in MeOH
Oven:	120 °C to 290°C at 15 °C min <sup>-1</sup> , hold 2min
Transfer Line:	280 °C

#### **Mass Spectrometry Conditions**

230 °C.
150 °C.
7 min.
SIM.
1 - 7.4 min to 7.8 min / 2 - 7.8 min to 8.1 min /
3 - 8.1 min to 8.4 min / 4 - 8.4 min to 9.8 min /

5 - 9.8 min to 10.4 min / 6 - 10.4 min to 16.0 min

#### Table 1. SIM parameters

Scan function	Compound	Quant Ion	1 <sup>st</sup> Qual Ion	2 <sup>nd</sup> Qual Ion	Dwell / ms
1	Butalbarbital	196	181	25	100
1	Butabarbital	169	211	37	100
2	Amobarbital	169	225	29	100
3	Pentobarbital	169	225	33	100
4	Secobarbital	196	181	25	100
5	Hexobarbital	235	81	27	100
6	Phenobarbital	232	175	33	100

## Results

Figure 2 shows the mass chromatograms for all the extracted barbiturates spiked at 10 ng/mL. Figure 3 shows average analyte recoveries ranging from 103-108% for barbiturates spiked at 10 ng/mL, with RSDs < 10% (n=3).







Figure 3. Average analyte recoveries of a range of Barbiturate analytes at 10 ng/mL (n=3).

## **Ordering information**

Part number	Description	Quantity
820-0140-C	ISOLUTE SLE + 1 mL sample	30

#### NORTH AMERICA

Main Office: +1 704 654 4900 Toll Free: +1 800 446 4752 Fax: +1 704 654 4917 Order Tel: +1 704 654 4900 press (4) at the auto attendant Order Fax: +1 434 296 8217 ordermailbox@biotage.com 1-pointsupport@biotage.com

#### EUROPE

Main Office: +46 18 56 5900 Fax: +46 18 59 1922 Order Tel: +46 18 56 57 10 Order Fax: +46 18 56 57 05 order@eu.biotage.com

#### Japan

Tel: +81 3 5627 3123 Fax:+81 3 5627 3121 jp\_order@biotage.com

#### Distributors

To locate a distributor please visit our Web site at www.biotage.com

 $\odot$  2011. All rights reserved. All brand and product names are trademarks or registered trademarks of their respective companies. The information contained in this document is subject to change without notice.

