

For Research Use Only. NOT for Use in Diagnostic Procedures.

Extraction of THC, Hydroxy-THC and Carboxy-THC from Whole Blood Using ISOLUTE® SLE+ Prior to GC/MS Analysis




Figure 1. Structure of Δ^9 -THC (tetrahydrocannabinol)

This application note describes the extraction of Δ^9 -THC, 11-hydroxy- Δ^9 -THC and 11-nor-9-carboxy-THC from whole blood matrix, prior to GC/MS analysis

Introduction

ISOLUTE® SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation.

This application note describes an effective and efficient ISOLUTE SLE+ protocol optimized for 1 mL capacity ISOLUTE SLE+ columns. The simple sample preparation procedure delivers clean extracts and analyte recoveries greater than 74% with RSDs lower than 9% for all analytes.

Analytes

THC, THC-OH, THC-COOH with THC-D3, THC-OH-D3 and THC-COOH-D3 as internal standards

Sample Preparation Procedure

- Sample Pre-treatment:** To 1.2 mL of whole blood, add 0.4 mL of 0.1% formic acid (aq), and mix thoroughly.
- Format:** **ISOLUTE® SLE+ 1 mL Sample volume columns, part number 820-0140-C**
- Sample loading:** Load 800 μ L of the pre-treated whole blood onto the column and apply a pulse of vacuum or positive pressure (3–5 seconds) to initiate flow. Allow the sample to absorb for 5 minutes.
- Analyte extraction:** Apply MTBE (3 mL) and allow to flow under gravity for 5 minutes. Apply Hexane (3 mL) and allow to flow for another 5 minutes under gravity. Apply vacuum or positive pressure (5–10 seconds) to pull through any remaining extraction solvent.
- Post elution and reconstitution:** Dry the extract in a stream of air or nitrogen using a SPE Dry (40 °C, 20 to 40 L/min) or TurboVap® (1.0 bar at 40 °C for 40 mins).
- Upon dryness, reconstitute with 40 μ L ethyl acetate and 20 μ L BSTFA:TMCS 99:1 and vortex for 20 seconds. Transfer to a high recovery glass vial. Place in a heating block set to 70 °C, for 25 minutes. Remove vial from the block and allow cooling.

GC Conditions

Instrument:	Agilent 7890A with QuickSwap
Column:	Agilent J&W DB-5, 30 m x 0.25 mm ID x 0.25 µm
Carrier:	Helium 1.2 mL/min (constant flow)
Inlet:	250 °C, Splitless, purge flow: 50 mL/min at 1.0 min
Injection:	2 µL
Wash solvents:	Acetone and Ethyl acetate
Oven:	Initial temperature 60 °C Ramp 25 °C/min to 350 °C, hold for 0.4 minutes
Post run:	Backflush for 2.4 minutes (3 void volumes)
Transfer Line:	280 °C

MS Conditions

Instrument:	Agilent 5975C
Source:	230 °C
Quadrupole:	150 °C
MSD mode:	SIM

SIM Parameters

Table 1. Ions acquired in the Selected Ion Monitoring (SIM) mode

SIM Group	Analyte	Target (Quant) Ion	1st Qual Ion	2nd Qual Ion
1	THC-D3	374	306	315
1	THC	317	343	386
2	THC-OH-D3	374		
2	THC-OH	371		
3	THC-COOH-D3	374	491	
3	THC-COOH	371	488	

Results

The optimized SLE+ protocol demonstrated analyte recoveries ranging from 72-83% as shown in **Figure 2**. RSD's were below 10% for all analytes.




Figure 2. Typical extraction % recoveries (n=7) using the ISOLUTE® SLE+ protocol.




Figure 3. Calibration curves for extracted levels of spiked whole blood using 1 mL capacity ISOLUTE® SLE+ columns from 1 ng/mL to 150 ng/mL showing r^2 values of 0.999 or greater.




Figure 4. GC/MS chromatography of whole blood. All analytes spiked at 1 ng/mL.

Table 2. Lower Limits of Quantitation (LLOQ) using this ISOLUTE® SLE+ procedure

Analyte	Lower Limit Of Quantitation
THC	1 ng/mL
THC-OH	3 ng/mL
THC-COOH	3 ng/mL

Additional Information

1. All solvents were HPLC grade.
2. 0.1% formic acid in water was prepared by adding 50 µL concentrated formic acid to 49,950 µL HPLC grade water.
3. Columns are slightly underloaded (800 µL of pre-treated sample on a 1 mL total capacity column) to ensure clean extracts from whole blood.

Ordering Information

Part Number	Description	Quantity
820-0140-C	ISOLUTE® SLE+ 1 mL Sample Volume Column	30
820-0140-CG	ISOLUTE® SLE+ 1 mL Sample Volume Column (tablets) for use on Biotage® Extrahera	30
820-0140-C-1000	ISOLUTE® SLE+ 1 mL Sample Volume Column (Bulk pack)	1000
820-0140-CG-1000	ISOLUTE® SLE+ 1 mL Sample Volume Column (tablets) (Bulk Pack) for use on Biotage® Extrahera™	1000
PPM-48	Biotage® PRESSURE+ 48 Positive Pressure Manifold for Columns	1
SD-9600-DHS-EU	Biotage® SPE Dry Sample Concentrator System 220/240 V	1
SD-9600-DHS-NA	Biotage® SPE Dry Sample Concentrator System 100/120 V	1
C103198	TurboVap® LV, Evaporator 100/120V	1
C103199	TurboVap® LV, Evaporator 220/240V	1

EUROPE

Main Office: +46 18 565900
 Toll Free: +800 18 565710
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
 order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: + 46 18 56 57 11
 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 704 654 4900
 Order Fax: +1 434 296 8217
 ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 Outside US: +1 704 654 4900
 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
 jp_order@biotage.com
 jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
 Fax: +86 21 2898 6153
 cn_order@biotage.com
 cn-1-pointsupport@biotage.com

To locate a distributor,
 please visit our website
 www.biotage.com

Part Number: AN840.V.1

© 2015 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe. **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**