For Research Use Only. NOT for Use in Diagnostic Procedures.

Extraction of Veterinary Growth Promoters from Animal Urine Using ISOLUTE[®] HYDRO DME+ Prior to UPLC-MS/MS Analysis



Figure 1. Representative veterinary growth promoter chemical structures.

This application note describes the extraction of 18 growth promoters (steroids, resorcylic acid lactones and stilbenes) using ISOLUTE® HYDRO DME+ prior to LC/MS analysis.

ISOLUTE HYDRO DME+ dual mode extraction columns and plates provide extremely efficient removal of matrix components and hydrolysis enzymes from urine samples, using a simple flowthrough workflow. This application note provides an effective and efficient ISOLUTE HYDRO DME+ procedure using column and 96-well plate formats. This simple sample preparation delivers clean extracts and recoveries greater than 80% for most analytes. RSDs are typically less than 5% and lower than 10% for all analytes. LLOQ are typically less than the estimated linear range minima of 0.4 ng mL⁻¹, using a total load volume of 75 µL of urine.

Analytes

β-trenbolone, α-trenbolone, β-boldenone, α-boldenone,
β-nortestosterone, α-nortestosterone, methylboldenone,
methyltestosterone, 16-β-hydroxystanozolol, taleranol, zeranol,
β-zearalenol, α zearalenol, zearalenone, zearalanone, diethylstilbestrol, dienestrol, hexestrol

Internal Standards

N/A

Sample Preparation Procedure

Format

ISOLUTE[®] HYDRO DME+ 400 mg columns, part number 970-0040-BZ

ISOLUTE° HYDRO DME+ 400 mg plate, part number 970-0400-PZ01

Matrix Preparation

Hydrolyze 75 μ L animal urine as follows: dilute 1:1 (v./v) with 200 mM ammonium acetate pH 5.2 containing 10 μ L per mL of glucuronidase H-2 ex. *Helix pomatia*. Incubate at 60 °C for 2 hours. Cool to room temperature.

Sample Loading

Load 150 μ L of hydrolyzed sample directly to ISOLUTE® HYDRO DME+ column/well. Add 900 μ L of acetonitrile and mix with 5x aspirate/dispense steps.

Analyte Extraction

Using a Biotage[®] Pressure+ 48 or 96 Positive Pressure Manifold, apply approximately 5 psi pressure. Collect the extract.

Post Elution and Reconstitution

Evaporate the extract at 40 °C and reconstitute in 0.1% formic acid in 50% aq MeOH (200 $\mu L).$



UHPLC Conditions

Instrument

Shimadzu Nexera UHPLC

Column

ACE UltraCore Super C18 2.5 µm 100 x 2.1 mm (CORE-25A-1002U), VWR International (Lutterworth, UK)

Mobile Phase

A: 0.5 mM ammonium acetate, 0.01% acetic acid, 30% MeOH (aq)

B: 0.5 mM ammonium acetate, 0.01% acetic acid, 95% MeOH (aq)

Flow Rate

0.5 mL min

Injection Volume

10 µL

Column Temperature

45 °C

MS Conditions

Instrument

Sciex Triplequad 5500 operating in dual polarity ESI mode

Source Temperature (TEM)

600 °C

Curtain Gas (CUR)

35

Source Gas 1 (GS1)

60

Source Gas 1 (GS2)

50

Table 3. MS parameters for target analytes.

Table 1. Gradient and divert valve settings Time (min) %A %B **Divert Valve** 0.00 10.0 90.0 7.00 NA NA 8 00 75 0 25 0

0.00	/ 5.0	23.0	
12.00	62.5	37.5	
14.00	NA	NA	waste
14.50	50.0	50.0	
16.50	10.0	90.0	
18.00	10.0	90.0	
20.00	90.0	10.0	
24.50	90.0	10.0	

waste

MS

Table 2. Typical retention times.

Analyte	Retention Time (min)
taleranol	7.5
β-trenbolone	8.1
β-zearalenol	8.4
β-boldenone	9.1
a-trenbolone	9.4
β-nortestosterone	9.6
zeranol	10.5
methylboldenone	10.6
a-zearalenol	11.2
zearalanone	11.5
a-boldenone	11.7
diethylstilbestrol	11.9
zearalenone	12.0
a-nortestosterone	12.0
16-β-hydroxystanozolol	12.1
dienestrol	12.6
methyltestosterone	12.9
hexestrol	13.1

Analyte	Transition (Da)	IS (V)	DP (V)	EP (V)	CE (V)	CXP (V)
taleranol	321.2 > 277.1	-3000	-112	-4	-30	-43
β-trenbolone	271.0 > 199.0	+4500	48	5	30	22
β-zearalenol	319.3 > 159.9	-3000	-85	-4	-38	-18
β-boldenone	287.2 > 135.0	+4500	73	9	21	18
a-trenbolone	271.0 > 199.0	+4500	30	4.5	32	22
β-nortestosterone	275.2 > 108.8	+4500	49	5.5	32	14
zeranol	321.2 > 277.0	-3000	-100	-4	-30	-38
methylboldenone	301.1 > 120.8	+4500	85	9.5	30	16.5
a-zearalenol	319.2 > 275.0	-3000	-180	-4	-28	-44
zearalanone	319.3 > 275.0	-3000	-90	-7	-28	-48
a-boldenone	287.2 > 120.9	+4500	62	7.5	31	18
diethylstilbestrol	267.1 > 237.0	-3000	-210	-8	-37	-48
zearalenone	317.2 > 131.0	-3000	-110	-4	-32	-17
a-nortestosterone	275.0 > 108.8	+4500	20	4.5	33	14
16-β-hydroxystanozolol	345.4 > 94.90	+4500	70	5	51	13
dienestrol	265.0 > 92.80	-3000	-130	-4	-31	-12
methyltestosterone	303.2 > 97.00	+4500	82	8	31	17
hexestrol	269.1 > 134.0	-3000	-74	-4	-20	-17



Results

Recovery

Extraction recoveries were determined at 1ng mL⁻¹, equivalent to 75 pg when extracting 150 μ L of pre treated hydrolyzed urine from four different species: bovine, equine, ovine and porcine.



Figure 2. Representative analyte recoveries using an optimized ISOLUTE® HYDRO DME+ protocol (75 µL of urine extracted).



Figure 3. Representative +ESI target analyte chromatography (1 ng mL $^{-1}$ spike).



Figure 4. Representative -ESI target analyte chromatography (1 ng mL 1 spike).



Linearity

Extraction linearity was determined between 0.4 and 10 ng mL⁻¹ from a mixed stock spiked into bovine urine and serially diluting in bovine urine. Each calibration level was extracted in four replicates. Figure 4 demonstrates representative calibration curves. Table 4 details linearity performance and associated LOQ for each analyte. Good linearity was observed for all analytes typically delivering r² values greater than 0.999.

Table 4. Analyte	calibration	curve r ²	and	LOQ	performance
------------------	-------------	----------------------	-----	-----	-------------

Analyte	Coefficient (r²)	Linear range, ng mL ⁻¹	LOQ, ng mL ⁻¹	% accuracy (RSD) at LOQ
taleranol	0.9975	0.4 to 10	0.4	95.8 (10.3)
β-trenbolone	0.9995	0.4 to 10	0.4	99.1 (3.7)
β-zearalenol	0.9970	0.4 to 10	0.4	97.1 (10.1)
β-boldenone	0.9994	0.4 to 10	0.4	98.1 (6.7)
a-trenbolone	0.9981	0.4 to 10	0.4	101 (5.8)
β-nortestosterone	0.9996	0.4 to 10	0.4	98.3 (3.5)
zeranol	0.9990	0.4 to 10	0.4	94.9 (2.7)
methylboldenone	0.9997	0.4 to 10	0.4	98.4 (4.9)
a-zearalenol	0.9982	0.4 to 10	0.4	99.8 (12.5)
zearalanone	0.9967	0.4 to 2.0	0.4	103 (6.5)
a-boldenone	0.9995	0.4 to 10	0.4	97.7 (1.1)
diethylstilbestrol	0.9969	0.4 to 10	0.4	102 (8.1)
zearalenone	0.9994	0.4 to 10	0.4	96.3 (7.4)
a-nortestosterone	0.9970	0.4 to 10	0.4	99.0 (6.0)
16-β-hydroxystanozolol	0.9984	0.4 to 10	0.4	95.9 (3.4)
dienestrol	0.9996	0.4 to 10	0.4	94.8 (5.9)
methyltestosterone	0.9996	0.4 to 10	0.4	99.8 (2.9)
hexestrol	0.9993	0.4 to 10	0.4	94.2 (4.0)



Figure 5. Representative Calibration curves (0.4 to 10 ng mL⁻¹ extracted from bovine urine) a) β -zearalenol, b) β -boldenone, c) diethylstilbestrol, d) methyltestosterone.



Chemicals and Reagents

Water (18.2 M Ω .cm) was drawn fresh daily from a Millipore Direct-Q5 water purifier (Watford, UK).

Organic solvents (methanol and acetonitrile) were HPLC or LC-MS grade and purchased from Honeywell Research Chemicals (Bucharest, Romania).

Acetic acid (LC-MS) was purchased from Fisher Scientific UK (Loughborough, UK).

Ammonium acetate (reagent and LC-MS grade), formic acid (LC-MS grade) and ß-glucuronidase (Type HP-2, aqueous solution ≥100,000 units/mL) were purchased from Sigma-Aldrich Company Ltd. (Gillingham, UK).

Analytical standards were purchased from Sigma-Aldrich Company Ltd. (Gillingham, UK) or Toronto Research Chemicals (Toronto, Canada)

200 mM ammonium acetate pH 5.2 (aq) was prepared by dissolving 15.42 g of reagent grade ammonium acetate in 1 L deionized water, the pH was adjusted using concentrated formic acid.

Buffered ß-glucuronidase hydrolysis solution was prepared fresh daily by diluting 150 μ L ß-glucuronidase solution in 14.85 mL 200 mM ammonium acetate pH 5.2 (aq) and vortexmixing thoroughly.

50% methanol (aq), 0.1% formic acid (reconstitution solution) was prepared by adding 50 mL methanol to 50 mL deionised water, followed by 100 μ L concentrated LC-MS grade formic acid.

Mobile phase A (0.5 mM ammonium acetate (aq), 0.01 % acetic acid) was prepared by dissolving 38.5 mg LC-MS grade ammonium acetate in 500 mL deionized water, adding 100 μ L of concentrated acetic acid and making up to 1 L with deionized water.

Mobile phase B (0.5 mM ammonium acetate (methanol), 0.01 % acetic acid) was prepared by dissolving 38.5 mg LC-MS grade ammonium acetate in 500 mL LC-MS grade methanol, adding 100 μ L of concentrated acetic acid and making up to 1 L with LC-MS grade methanol.

Additional information

Animal urine matrix samples were kindly donated by the State Laboratory, Celbridge, Ireland.

Ordering Information

Part Number	Description	Quantity
970-0400-PZ01	ISOLUTE [®] HYDRO DME+ 400 mg Plate	1
970-0040-BZ	ISOLUTE® HYDRO DME+ 400 mg/3 mL Columns	50
PPM-96	Biotage [®] PRESSURE+ 96 Positive Pressure Manifold	1
PPM-48	Biotage PRESSURE+ 48 Positive Pressure Manifold	1
45000	TurboVap® LV	1
SD-9600-DHS-NA	SPE Dry 96 Sample Concentrator System, 100/120V	1
SD-9600-DHS-EU	SPE Dry 96 Sample Concentrator System, 220/240V	1

EUROPE

Main Office: +46 18 565900 Toll Free: +800 18 565710 Fax: +46 18 591922 Order Tel: +46 18 565710 Order Fax: +46 18 565705 order@biotage.com Support Tel: +46 18 56 59 11 Support Fax: +46 18 56 57 11 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA Main Office: +1 704 654 4900

Toll Free: +1 800 446 4752 Fax: +1 704 654 4917 Order Tel: +1 704 654 4900 Order Fax: +1 434 296 8217 ordermailbox@biotage.com Support Tel: +1 800 446 4752 Outside US: +1 704 654 4900 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123 Fax: +81 3 5627 3121 jp_order@biotage.com jp-1-pointsupport@biotage.com

CHINA Tel: +86 21 68162810 Fax: +86 21 68162829 cn_order@biotage.com cn-1-pointsupport@biotage.com

KOREA Tel: +82 31 706 8500 Fax: +82 31 706 8510

Fax: +82 31 706 8510 korea_info@biotage.com kr-1-pointsupport@biotage.com

INDIA Tel: +91 22 4005 3712 india@biotage.com

Distributors in other regions are listed on www.biotage.com

Ø Biotage

Literature Number: AN925.V.1

© 2019 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe. **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**