


For Research Use Only. NOT for Use in Diagnostic Procedures.

# Improving Analytical Sensitivity and Sample Throughput Using Biotage® Mikro Solid Phase Extraction Microelution Plates

## Introduction


Solid phase extraction is a powerful extraction technique which delivers high analyte recoveries and clean extracts. Microelution SPE formats utilize very small amounts of extraction sorbents, suitable for processing low volumes of sample, and can deliver excellent analyte recovery in very low elution volumes.

Biotage® Mikro solid phase extraction microelution plates have a very low dead volume design, with relatively tall, narrow beds. These dimensions mean that the Mikro plates have high analyte capacity, but analyte recovery is achieved using much lower elution volumes than traditional SPE formats.



**Figure 1.** Comparison of SPE bed dimensions 10 mg SPE plate vs Biotage® Mikro Plate (not to scale).

The combination of high capacity and low elution volume leads to an increased concentration of analyte in the elution solvent.



**Figure 2.** Schematic of workflow using 10 mg SPE plate vs Biotage® Mikro plate.

## Experimental

In this example, we compared the use of Biotage® Mikro ABN plates with traditional format EVOLUTE® EXPRESS ABN 10 mg 96-well plates for the extraction of steroids from human serum. Serum was spiked at a concentration of 1 ng/mL. The same extraction sorbent (EVOLUTE ABN) was packed into both formats, and the same methodology was used in each case. Solvent volumes were reduced for the Mikro format (see table 1 for the SPE methods used). Both plates were processed using the Biotage® PRESSURE+ 96 Positive Pressure Manifold.

After extraction, the plate was sealed with a sealing mat until analysis. 5 µL of diluted extract was analysed by LC-MS/MS (see Biotage application note AN939 for full analytical details), and analyte sensitivity was compared.


**Table 1.** SPE Methods.

Step	EVOLUTE® EXPRESS ABN 10 mg Plate	Biotage® Mikro ABN 2 mg Plate
<b>Condition</b>	Methanol (500 µL)	Methanol (100 µL)
<b>Equilibrate</b>	0.1% formic acid (500 µL)	0.1% formic acid (100 µL)
<b>Load</b>	400 µL diluted serum 200 µL serum diluted with 200 µL 1% formic acid	400 µL diluted serum 200 µL serum diluted with 200 µL 1% formic acid
<b>Wash 1</b>	Water (500 µL)	Water (100 µL)
<b>Wash 2</b>	H <sub>2</sub> O:MeOH (60:40, v/v, 500 µL)	H <sub>2</sub> O:MeOH (60:40, v/v, 100 µL)
<b>Elute</b>	Methanol (150 µL)	Methanol (30 µL)
<b>Dilute*</b>	Water (150 µL)	Water (30 µL)
<b>Inject</b>	5 µL	5 µL

\*Note after elution, the eluent was diluted (1:1, v/v) with water for compatibility with the analytical HPLC conditions.

## Results and Discussion

Using the methods described, very clean extracts were produced and excellent recovery with low RSDs was achieved using both SPE formats. Extracts produced using the Mikro format showed a significant gain in sensitivity compared to the traditional 10 mg format (see figure 3).



**Figure 3.** Extracted ion chromatogram for estradiol extracted from human serum using the methods outlined in table 1. Microelution SPE using Biotage® Mikro ABN plates resulted in a ~8 x increase in signal/noise.

## Conclusion

Biotage® Mikro ABN solid phase extraction microelution plates provide high analyte recoveries with excellent reproducibility. Compared to the traditional 10 mg SPE plate format, improved sensitivity and reduced LLOQs have been demonstrated. In order to achieve equivalent sensitivity using the 10 mg plate format, an evaporation and reconstitution step would be required, adding approximately 15 minutes to the overall sample preparation time and introducing the possibility of processing errors.

Low volume microelution can lead to immediate improvements in analyte sensitivity, without the need for an evaporation step.

## Ordering Information

Part Number	Description	Quantity
<b>600-0002-LVP</b>	Biotage® Mikro ABN Plate, 2 mg	1
<b>601-0002-LVP</b>	Biotage® Mikro CX Plate, 2 mg	1
<b>602-0002-LVP</b>	Biotage® Mikro WCX Plate, 2 mg	1
<b>603-0002-LVP</b>	Biotage® Mikro AX Plate, 2 mg	1
<b>604-0002-LVP</b>	Biotage® Mikro WAX Plate, 2 mg	1
<b>Collection Plate and Accessories</b>		
<b>121-5202</b>	Collection Plate, 1 mL Square	50
<b>121-5204</b>	Pierceable Sealing Mat	50
<b>Positive Pressure Manifold</b>		
<b>PPM-96</b>	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1

### EUROPE

Main Office: +46 18 565900  
Toll Free: +800 18 565710  
Fax: +46 18 591922  
Order Tel: +46 18 565710  
Order Fax: +46 18 565705  
order@biotage.com  
Support Tel: +46 18 56 59 11  
Support Fax: +46 18 56 57 11  
eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
Support Tel: +1 800 446 4752  
Outside US: +1 704 654 4900  
us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 68162810  
Fax: +86 21 68162829  
cn\_order@biotage.com  
cn-1-pointsupport@biotage.com

### KOREA

Tel: +82 31 706 8500  
Fax: +82 31 706 8510  
korea\_info@biotage.com  
kr-1-pointsupport@biotage.com

### INDIA

Tel: +91 22 4005 3712  
india@biotage.com

Distributors in other regions are listed on [www.biotage.com](http://www.biotage.com)

Literature Number: PPS619.V.1

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe. **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**