For Research Use Only. NOT for Use in Diagnostic Procedures.

Extraction of Fat-Soluble Vitamins from Human Serum Using ISOLUTE[®] SLE+ Prior to UHPLC/MS-MS Analycic



Figure 1. Structures of Retinol and Beta Carotene (Vitamin A), 25-OH-Vitamins D2 and D3 (Vitamin D), Alpha Tocopherol (Vitamin E) and Phylloquinone and Menaquinone-4 (Vitamin K).

This application note describes the extraction of a panel of fat-soluble vitamins (including those representing Vitamins A, D, E & K) from human serum using ISOLUTE[®] SLE+ Supported Liquid Extraction plates prior to LC/MS analysis.

The simple sample preparation procedure delivers clean extracts and analyte recoveries approximately or above 90% with RSDs lower than 10% for all analytes.

Analytes

Retinol, Beta Carotene, 25-OH Vitamin D2, 25-OH Vitamin D3, Alpha Tocopherol, Phylloquinone, Menaquinone-4.

Internal Standards

D₆ 25-OH Vitamin D3 was used as an internal standard for 25-OH Vitamin D2, 25 OH Vitamin D3 and Retinol. It is recommended that an additional internal standard is used for Vitamin K.

Table 1. Concentration ranges for fat-soluble vitamins.

Analyte name and Vitamin reference	Lower limit of Quantification (ng/mL)	Upper limit of Quantification (ng/mL)
Alpha Tocopherol (Vitamin E)	800	40000
Retinol and Beta Carotene (Vitamin A)	80	4000
25-OH Vitamins D2 & D3 (Vitamin D)	4	200
Phylloquinone and Menaquinone-4 (Vitamin K)	0.4	20



Sample Preparation Procedure

Format

ISOLUTE° SLE+ 400 μL supported liquid extraction plate, part number 820-0400-P01.

Sample Pre-treatment

Internal standard solution (10 μ L) and 10 μ L of either spiking solvent or calibration standard spiking solution (see 'Chemicals and Reagents' section for preparation details) were transferred to a 2 mL collection plate. 100 μ L of sample or blank matrix was added, capped, briefly mixed and then left to stand in the dark for 1 hour to equilibrate. The sample was then combined with 400 μ L of pre-treatment solvent, briefly mixed and then left to stand for a further 5 minutes.



Figure 2. Demonstrating the different preparation procedure for samples and Calibration Standards.

Sample Loading

Approximately 500 μ L pre-treated serum (or as much of the sample as possible) was transferred to an ISOLUTE^{*} SLE+ plate. The sample was vigorously drawn up and down into the pipette tip with additional air (by setting a larger volume than the sample e.g. 650 μ L) a number of times immediately prior to transfer to create a temporary suspension. If necessary, a low positive pressure was used to push the sample into the SLE material. The ISOLUTE SLE+ plate was then left to equilibrate for 5 minutes.

Elution

Analytes were eluted with 2 x 500 μ L heptane.

Post Elution and Reconstitution

The extract was dried in a stream of air or nitrogen using a Biotage[®] SPE Dry 96 at room temperature, 20 to 40 L/min.

Evaporated samples were reconstituted with propan-2-ol (IPA, 150 $\mu L)$ and mixed thoroughly.

UPLC Conditions

Instrument

Waters Acquity UPLC

Column

Restek Raptor Biphenyl (100 mm x 2.1 mm, 2.7 $\mu\text{m})$ with a Restek EXP holder and guard.

Mobile Phase

A: 5 mM ammonium acetate 0.1% formic acid (v/v) in water

B: Methanol : Propan-2-ol (3:1 v/v) containing 5 mM ammonium acetate and 0.1% formic acid

Flow Rate

o.4 mL/min

Column Temperature

40 °C

Autosampler Temperature

10 °C

Injection Volume

10 μL (Partial Loop with Needle Overfill)

Table 2. UHPLC Gradient.

Time (min)	% A	%В
0	60	40
3	0	100
5.1	0	100
5.2	60	40
7	60	40



Biotage[®] SPE Dry Sample Concentrator System.



MS Conditions

Instrument Waters Quattro Premier XE

Desolvation Gas Flow 1200 L/hr

Cone Gas Flow

Source Temp

Desolvation Temp 450 °C

Capillary Voltage 4 kV

Extractor Voltage

Results

Extraction recoveries were first measured using a manual processing method (using a Biotage[®] PRESSURE+ 96 manifold). The method was then transferred to a Biotage[®] Extrahera[®] for automated processing. The Extrahera[®] recoveries were slightly lower in line with the slower mixing of this compared to the manual method. Extraction recoveries (manual and automated methods), and associated RSDs are shown in table 4.

Table 3. MS conditions and retention times for target analytes.

Analytes	MRM Transition	Collision Energy, V	Cone, V	Period
Alpha Tocopherol (E)	433.3 > 165.9	22	25	2
Retinol (A)	269.3 > 92.9	20	18	1
Beta Carotene (A)	536.3 > 444.4	15	30	2
25-OH Vitamin D2	395.5 > 269.5	30	30	1
25-OH Vitamin D3	383.5 > 257.5	17	30	1
Phylloquinone (K1)	445.3 > 186.9	20	22	2
Menaquinone-4 (K2)	451.4 > 187.0	23	30	2
D ₆ 25-OH Vitamin D3 (IS)	389.6 > 263.5	16	30	1

All analytes were measured in positive mode using Electrospray ionization.

Note: The transition for Alpha Tocopherol was significantly different from the optimum settings. Due to the high MS sensitivity of the analyte and the high concentrations expected this was intentionally de-tuned on the instrument.

Table 4. Analyte calibration curve r² and LOQ performance.

Analytes	Manual		Biotage [®] Extrahera™		
	Recovery	% RSD	Recovery	% RSD	
Retinol (A)	94.9	3.3	75.9	5.3	
Beta Carotene (A)	89.0	4.1	71.3	10.5	
25 OH Vitamin D2 (D)	101.0	5.1	79.2	5.4	
25 H Vitamin D3 (D)	95.6	4.6	81.9	5.5	
Alpha Tocopherol (E)	99.1	4.6	84.0	4.6	
Phylloquinone (K)	95.7	10.5	71.0	6.3	
Menaquinone-4 (K)	95.7	9.6	73.4	4.5	

 D_6 25-OH Vitamin D3 was used as an internal standard for 25-OH Vitamin D2, 25-OH Vitamin D3 and Retinol.



Biotage[®] PRESSURE+ 96 Positive Pressure Manifold.









Figure 4. Representative chromatography for stripped serum spiked at a midcalibration range (6 ng/mL Vitamin K, 60 ng/mL Vitamin D, 1.2 µg/mL Vitamin A and 12 µg/mL Vitamin E.



Discussion and Conclusion

This method provides high, reproducible recoveries of a range of fat-soluble vitamins in human serum, in clinically appropriate concentration ranges.

Due to the extremely non-polar (hydrophobic) nature of the analytes, and the wide difference in biological concentration range, some non-standard modifications to the standard supported liquid extraction process were adopted in this application. See the 'Additional Method Notes' section for a comprehensive description of the steps taken to ensure successful extraction of these analytes.

Chemicals and Reagents

Stock and Sub Stock Solvent Where solids were provided these were diluted in a stock solution solvent of MTBE + 1 mg/mL BHT.

Stock solutions were prepared at 1 mg/mL (Retinol, Phylloquinone, Menaquinone). Due to the challenges of precisely weighing an oil Alpha Tocopherol was prepared in stock solution solvent at known concentrations between 1 and 5 mg/mL. Due to limited solubility beta carotene was prepared in stock solution solvent at a concentration of 100 μ g/mL. 25-OH Vitamin D2 and D3 were purchased as solutions. All stock solutions were stored protected from light at approximately -20 °C.

Spiking solvent

Spiking solvent was prepared by combining BHT in propan-2-ol (IPA) at a level of 0.1% w/v or 1 mg/mL. As an example: 100 mL of propan-2-ol would be added 100 mg of BHT.

Internal standard solution

Internal standard (D⁶ 25-OH Vitamin D3) was diluted in spiking solvent to a concentration of 1 μ g/mL. A 10 μ L aliquot of this is equivalent to 100 μ L of sample containing internal standard at a level of 100 ng/mL.



Combined Fat-Soluble Vitamin Spiking Solution

A spiking solvent was prepared by combining fat-soluble vitamin solutions and diluting with spiking solvent such that the following concentrations were met: Vitamin K1 and K2D4 = 200 ng/mL, 25-OH Vitamin D2 and D3 = 2 μ g/mL, Retinol and Beta Carotene = 40 μ g/mL and Alpha Tocopherol = 400 μ g/mL. It is recommended that at least 0.5 mL of this solution is prepared on a daily basis. The "Additional Information" section contains an example spiking procedure to reach the required concentrations. It is recommended that this solution is prepared daily.

Calibration standards

Calibration standard spiking solutions were prepared from the combined fat-soluble vitamin spiking solution. The "Additional Information" section contains an example spiking procedure to reach the required concentrations. It is recommended that this solution is prepared daily.

Pretreatment solvent

BHT was combined with a solution of IPA/heptane (1:3, v/v) at a level of 1 mg/mL.

Elution solvent

Heptane was used as the SLE elution solvent. Hexane is an acceptable alternative, but extract cleanliness may be slightly compromised (slightly higher levels of co-extracted phospholipids may be observed).

Reconstitution Solvent

Propan-2-ol was used as the reconstitution solvent.

Other Chemicals and Reagents

- » Method development was performed using vitamin stripped serum was purchased from Golden West.
- Methanol (LC-MS grade), propan-2-ol (isopropanol) (99.9%), MTBE (99%) ethyl acetate and formic acid (98%) were purchased from Honeywell Research Chemicals (Bucharest, Romania).
- » All analyte standards and deuterated internal standards were purchased from Sigma- Aldrich Company Ltd. (Gillingham, UK).
- Water used was 18.2 MOhm-cm, drawn from a Direct-Q5 water purifier.
- » Mobile phase A was prepared by accurately weighing approximately 385.4 mg of ammonium acetate. The ammonium acetate (385.4 mg assumed) was then combined with 1 L of water and 1 mL of formic acid. The solution was replaced after 48 hours.
- Mobile phase B was prepared by accurately weighing approximately 385.4 mg of ammonium acetate. The ammonium acetate (385.4 mg assumed) was then combined with 750 mL of methanol, 250 mL of propan-2-ol and 1 mL of formic acid. The solution was replaced after 48 hours.
- Mobile Phase B was partially prepared with propan-2-ol to give improved chromatographic retention times of the later eluting analytes, particularly Beta-Carotene.
- Due to analyte instability the preservative BHT is included in all stock and pre-treatment solvents. BHT is also known as 2,6 Di-tert-butyl-4-methylphenol (Sigma Aldrich (B1378)).



Additional Method Notes

- In vivo most fat-soluble vitamins are highly protein bound » so the serum samples must be treated to significantly reduce these interactions prior to any extraction.
- » Any sample incubations were performed in the dark. If using the Biotage[®] Extrahera[®] the internal lights were switched off.
- » By adding heptane as a pre-treatment solvent, the solubility of protein freed fat-soluble vitamins in the pre-treated sample was increased. Rapidly aspirating and dispensing the sample in the tip prior to analysis formed a temporary milky colored suspension. Suspension formation was aided by the presence of propan-2-ol which is separately soluble in both serum and heptane.



Figure 5. Images representing sample consistency of sample alone, during sample incubation, after addition of pretreatment solvent and immediately after mixing.

- A volume of 500 µL could be loaded on to the » ISOLUTE® SLE 400 because the sample only contained approximately 200 µL of aqueous solvent.
- » Due to the consistency of the pre-treated serum it may not flow into the SLE bed under gravity and may require application of positive pressure or vacuum for loading.
- » Increased vitamin losses were seen when the extracted samples were evaporated at 40 °C and so this step is performed at room temperature.
- » Propan-2-ol was used as the reconstitution solvent in order to keep the beta carotene in solution over the period of the LC analysis.
- » The method detailed here involved the use of one internal standard that was used for the measurement of 25-OH Vitamin D2, 25-OH Vitamin D3 and Retinol. Due to the suppression levels seen with Vitamin K measurement, the use of an isotope or structurally similar compound as an additional internal standard to Vitamin K1 and K2D4 is strongly recommended.
- » The calibration line for Alpha Tocopherol will show an element of non-linearity due to the high levels of analyte being injected. A similar non-linear plot may also be present for Beta Carotene.
- The preparation method of quality control samples should » be carefully considered due to the possible binding of some fat-soluble vitamins to polypropylene in a highly aqueous environment. Although this isn't an issue for the analytical procedure described here, as the analytes would come back into solution once heptane had been mixed with it, it could be an issue for externally prepared QCs.

- **>>** All calibration ranges were set with the aim of quantifying the majority of normal samples. With the Vitamin K analytes these ranges could not be confirmed and so the lowest range that could be confidently measured was used. Due to the small Vitamin K peaks generated it is recommended that if the experiment is performed in an area where manual integration is not permitted, the extracts are analyzed on a more sensitive MS than the one used in this experiment.
- » Vitamin K1 and K2D4 were the only compounds that showed a significant level of suppression. This was shown to decrease significantly when less serum was extracted. If running the method on a more sensitive LC-MS than is documented here it is recommended that a reduction of sample volume from 100 to 50 or 20 μ L is considered. Note if this is performed the concentration of analytes in sub stock 2 will need to be reduced accordingly.
- » To optimize MS performance the analyte detection was separated into two periods. Vitamin Ds (D3, D2 and D6-D3 (IS)) and Retinol were measured in period 1 from 0 to 3.5 minutes. All other components were measured between 3.5 and 6 minutes. »
- For increased sensitivity:
 - » Decrease reconstitution solvent volume below 150 µL
 - Consider using an LC-MS instrument >> with greater sensitivity
 - >> Due to the variability of stock solution concentrations and expected levels of each component in serum, the preparation of appropriate calibration range standards can be challenging. Below is an example of the preparation of calibration standards in a quantity appropriate for up to n=4 analysis which can be used as a guide: Stock solution concentrations assumed.
 - Retinol = 1 mg/mL»
 - Beta carotene = $100 \,\mu g/mL$ »
 - 25-OH Vitamin D2 = 50 μ g/mL »
 - 25-OH Vitamin D3 = 100 μ g/mL »
 - » Alpha Tocopherol = 2 mg/mL
 - Phylloquinone = 1 mg/mL»
 - Menaquinone-4 = 1 mg/mL »
 - » Substock 1. Combine 10 µL of Phylloquinone (K1) with 10 µL Menaquinone-4 (K2D4) and then dilute this mixture to a total volume of 1 mL with the addition of stock and sub stock solvent (MTBE + 0.1% w/v BHT). This solution contains Phylloquinone and Menaquinone at concentrations of 10 μ g/mL each.
 - Substock 2. Combine 10 µL of substock 1, 10 µL of >> 25-OH Vitamin D3, 20 µL of retinol, 20 µL of 25-OH Vitamin D2, 100 µL of Alpha Tocopherol and 200 µL of Beta Carotene. Dilute this mixture to a total volume of 0.5 mL with precipitation solvent (IPA + 0.1% BHT).
 - >> This solution contains Alpha Tocopherol at $400 \ \mu g/mL$, Retinol and Beta Carotene at 40 µg/mL, 25-OH Vitamin D2 and D3 at 2 μ g/mL and Vitamin K1 and K2 at 0.2 μ g/mL.

Calibration standards were prepared by combining substock 2 with spiking solvent (see Table 5).



Table 5. Spiking regimen of calibration standard spiking solutions.

Standard ID	Volume of Substock 2 (µL)	Volume of Spiking Solvent (µL)
1	10	490
2	10	240
3	20	230
4	40	210
5	30	70
6	60	40
7	100	0

Ordering Information

Part Number	Description	Quantity
820-0400-P01	ISOLUTE $^{\circ}$ SLE+ 400 μ L Supported Liquid Extraction Plate	1
PPM-96	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1
SD-9600-DHS	Biotage [®] SPE Dry Sample Concentrator system	1
121-5203	Collection Plate, 2 mL Square	50
121-5204	Pierceable Sealing Mat	50



Appendix Biotage® Extrahera™ Settings

The method described in this application note was automated on the Biotage[®] Extrahera[¬] using ISOLUTE[®] SLE+ 400 μ L capacity 96-well plates. This appendix contains the software settings required to configure Extrahera[¬] to run this method. Screenshots may or may not match those here depending upon the instrument software version.

Sample Name:	Fat Soluble Vitamins in Serum
Sample Plate/Rack:	2 mL 96 well FSV
Extraction Media:	ISOLUTE° SLE+ 400 µL 96







Settings "Sample" Tab

Sample Type: Starting Sample Volume (µL): Method Comment:

FSV Sample Mix 300

Prior to running on the Extrahera, combine 100 μ L serum, 10 μ L IS solution and 10 μ L Propan-2-ol or spiking solution then leave for 1hr to equilibrate. A larger than true starting sample volume has been set to maximise sample mixing prior to loading.

Pretreatment	
Steps	2
Solvent	Heptane : IPA (3:1)
Volume	400
Mix	5
Volume	700
Wait	5
Pause	No
Dispose tips	No



Fat Soluble V	itamin <mark>s in S</mark> e	erum	Sa	imple plate/rack 2 mL Sample Plat	te, 96 👻	Extraction media
Pretreatment	Sample	Pretreatment	Load	Elution (2)		
On .oad	Sample vo	lume (µL)		Air push time (s)	Advanced p settings	ressure
On	700			0	Edit.	
lution	Premix?	Number of	ftimes	Wait time (min)		
On	Yes	5	*	5		
	Pause afte	collect in r	position			
		No D (Wa.	👻			
			_			

Load	
Volume	700
Premix	Yes
No of times	5
Wait time	5
Pause	No
Collect in	D
Pressure	Advanced
Steps	2
0 bar	120
2 bar	30

Method name		Sample plate/rack	Extraction media
Fat Soluble Vit	amins in Serum	2 mL Sample Plate, 96 👻	ISOLUTE SLE+ 400 μL 9 👻
Pretreatment	Sample Pretreatment	Load Elution (2)	
On	Number of steps	Air push after last elution? Air push time	Dispose solvent tig
Load	2 -	No 0	No
On	1	(² column)	
Elution	Heptane	 Heptane 	-
On	Volume (µL) Collect în po	isition Volume (µL) Collect in positi	ion
	500 A	▼ 500 A	-
	Advanced p Wait time (min) Settings	verssure Wait time (min) Advanced pressettings	ssure
	0 Edit	t 0 Edit	
	Repeat (number Pause after	this Repeat (number Pause after th	lis
	of times) step?	No 1	No
Flution 1 A	dvanced settings	Flution 2 A	dvanced settings
Use advanced	Number of steps	Use advanced	Number of steps
Voc		pressure secong	
res	2	res	3 ×
Pressure (bar)	Positive pressure time (s)	Pressure (bar)	Positive pressure time (s)
1 0.0	180	0.5	180
Pressure (bar)	time (s)	Pressure (bar)	Positive pressure time (s)
2 2.0	30	2.0	60
L			Docitive accessure
Air Push?	Air push time (s)	Pressure (bar)	time (s)
N	0	5.0	10
		Air Push?	Air push time (s)
		No	0

Elution	Activated
Number of steps	2
Air push	No
Dispose tips	No
Air Push	No

Solvent 1 - Heptane	
Volume	500
Collect position	А
Advanced settings	Yes
Number of Adv. steps	2
Pressure 1	0
Time 1	180
Pressure 2	2
Time 2	30
Air push	No
Repeat	1
Pause	No

Solvent 2 - Heptane	
Volume	500
Collect position	A
Advanced settings	Yes
Number of Adv. steps	3
Pressure 1	0.5
Time 1	180
Pressure 2	2
Time 2	60
Pressure 3	5
Time 3	10
Air push	No
Repeat	1
Pause	No
Dispose tips	No



Solvent Properties

Solvent Description	
Heptane : Propan-2-ol	
Heptane	
	THE CALL
	TUTTE
	Solvent Description Heptane : Propan-2-ol Heptane

Solvent	1	2	3	4	5	6	7	8	9	10
Reservoir Type		Refillat	ole				No	on Refillable		
Capacity										
Aspiration flow rate	10	10								
Dispense flow rate	10	10								
Aspiration post dispense?	No	No								
Lower air gap flow rate	10	10								
Lower air gap volume	5	5								
Upper air gap flow rate	120	20								
Upper air gap volume	100	100								
Upper air gap dispense pause	200	200								
Conditioning?	Yes	Yes								
Frequency	1 st Asp.	only								
Cond. Times	4	4								
Cond. Flow rate	20	20								
Cond. volume	100	100								
Chlorinated	No	No								
Serial dispense	No	No								
Highly Volatile	No	No								

Sample	Air Gap
Sample name	Lower air ga
FSV Sample Mix	20
Sample description	Lower air ga
Fast mix, zero pause	5
Aspiration flow rate (mL/min)	Upper air ga
100	100
Dispense flow rate (mL/min)	Upper air ga
220	700
Aspirate post dispense?	Upper air ga
Yes	0
Aspirate post dispense flow rate (mL/min)	
20	
Aspirate post dispense volume (µL)	
100	
Territoria de la constante de la const	

20	
ower air gap v	volume (µL)
5	
pper air gap t	flow rate (mL/min)
100	
pper air gap v	volume (µL)
700	
pper air gap (dispense pause (ms)
0	

"Sample" Screen	
Sample name	Aqueous sample
Sample description	Fast mix, zero pause
Aspiration flow rate	100
Dispense flow rate	220
Aspirate post dispense?	Yes
Aspirate post dispense flow rate	20
Aspirate post dispense volume	100
Lower air gap flow rate	20
Lower air gap volume	5
Upper air gap flow rate	100
Upper air gap volume	700
Upper air gap dispense pause	0



Extraction Media	Pipetting Height
Name	Solvent dispensation height (mm)
ISOLUTE SLE+ 400 µL 96 Well	-125.0
Manufacturer	Sample dispensation height (mm)
Biotage	-126.5
Part number	Aspiration height (mm)
820-0400-P01	-139.5
Capacity volume (µL)	
0	Tune Pipetting Heights
Format	
96 Positions, Plate/Col 👻	
Comment	

"Extraction Media" Screen	
Name	ISOLUTE° SLE+ 400 µL 96 Well Plate
Manufacturer	Biotage
Part number	820-0400-P01
Capacity volume	0
Format	96
Comment	
Solvent dispensation height	-125
Sample dispensation height	-126.5
Aspiration height	-139.5

- <mark>1</mark> 61.5
Pretreatment dispensation height (mm)
-128.0
Tune Pinetting Heights
1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.

"Sample Plate/Rack" Screen	
Name	2 mL 96 well FSV
Capacity volume	1800
Format	96
Aspiration height	-161.5
Pretreatment dispensation height	-128



< Cancel	Edit Pipette Tip - 1000 µL Biotage tip	Save >
	Pipette Tip Name 1000 µL Biotage tip Manufacturer Biotage Part number 414141 1000 Length (mm) 95	

"Pipette tip" Screen	
Name	1000 µL Biotage Tip
Manufacturer	Biotage
Part number	414141
Capacity (µL)	1000
Length (mm)	95

EUROPE

Main Office: +46 18 565900 Toll Free: +800 18 565710 Fax: +46 18 591922 Order Tel: +46 18 565710 Order Fax: +46 18 565705 order@biotage.com Support Tel: +46 18 56 59 11 Support Fax: + 46 18 56 57 11 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900 Toll Free: +1 800 446 4752 Fax: +1 704 654 4917 Order Tel: +1 704 654 4900 Order Fax: +1 434 296 8217 ordermailbox@biotage.com Support Tel: +1 800 446 4752 Outside US: +1 704 654 4900 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123 Fax: +81 3 5627 3121 jp_order@biotage.com jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810 Fax: +86 21 68162829 cn_order@biotage.com cn-1-pointsupport@biotage.com

KOREA Tel: +82 31 706 8500 Fax: +82 31 706 8510 korea_info@biotage.com kr-1-pointsupport@biotage.com

INDIA Tel: +91 22 4005 3712 india@biotage.com

Distributors in other regions are listed on www.biotage.com



Literature Number: AN945.V.1

© 2020 Blotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explantion and to the owners' benefit, without intent to infringe. FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.