


Accelerating High Throughput Protein Purification Using PhyTip® Columns

A White Paper from Biotage

Author: David Westman



Introduction

Drug discovery and development is a multifaceted process which involves a long and exhaustive journey through basic research, discovery of the medicine, preclinical development tests, increasingly complicated clinical trials, and finally the regulatory approval. The initial drug discovery and development process includes target validation, high throughput screening of drug candidates, developing assays for therapeutic efficacy, and lead optimization. One of the key aspects in this phase is to become more efficient in finding the hits and shorten the time to identified drugs for process development and scaling up.

For biologics and advanced therapeutics, the purification steps for characterization and drug interaction studies can be a bottleneck. The process for adequate preparation often requires that sufficient quantities of material be scaled up and processed in a time-consuming manner using expensive chromatography equipment. Commonly used purification techniques are magnetic beads and filter plates. However, these technologies suffer from drawbacks such as expensive specific instrumentation for magnetic beads and complex automation protocols for the filter plates. Recent advancements in miniaturized high-throughput tools for purification, enrichment and desalting of proteins eliminate bottlenecks associated with traditional protein purification techniques.




Figure 1. The drug discovery and development process.

Dual Flow Chromatography

Dual flow chromatography (DFC) is the unique process of performing chromatography separations in a column bed loosely packed at the end of a pipette tip, retained by a minimally absorbing, proprietary mesh – a PhyTip® column. Dual flow chromatography works by having a flow through the PhyTip® column and pumping the sample back and forth across the packed resin bed inside the tip. As the sample mixture is cycled back and forth more than once with a controlled flow across the resin bed, the biomolecules in the sample is given enough time to interact with the available binding sites despite the slow transfer kinetics associated with large biomolecules. As a result, all the protein, viral vector or plasmid can be captured during the capture step, leading to very high recoveries at the end of the purification process. And if this is not the case, this is easily adjusted by simply applying more DFC cycles to increase the contact time between the sample and the binding sites on the resin. The same applies to the wash step and elution step (Figure 2) – the number of cycles is an easy way of fine-tuning the purification process for a specific application. And when each step of the protein purification protocol is completed this way and driven to equilibrium, the purification will perform the same every time. That means the results are highly reproducible with the same sample purity and yield every time.

The unique column design with a packed resin bed binds the entire sample to a small resin bed volume. Basic chromatographic principals control the separations, and the bi-directional flow drives interactions to high loading equilibrium regardless of kinetic rate constants. The purified sample can be eluted in small volumes, resulting in a highly concentrated sample. Dual flow chromatography is a gentle purification process producing high levels of biologically active protein.

Automated High Throughput Purification

Biotage' PhyTip® column technology has been developed for small scale high-throughput preparation of biomolecules such as antibodies, recombinant proteins, plasmids and viral vectors in order to facilitate the process of preparing hundreds to thousands of potential leads that are ready for cell-based assays without the need for scale-up. These unique columns are designed to operate on all major liquid handling robot platforms (Hamilton, Tecan, Beckman and more) either on 96-at-a-time platforms or on 8-12 channel robotic platforms, see Figure 3 below. The PhyTip® columns are ready to use directly out of the box. The run time for a full purification protocol of 96 samples can be as short as 15 minutes depending on sample volume. Methods are easy to develop, and validated protocols are available for many of these platforms.




Figure 2. Principles of dual flow chromatography.





Figure 3. PhyTip® columns are designed to operate on all major liquid handling robots.

Enabling High-Throughput Screening


When automating library screens, researchers have several goals they need to accomplish. They have to express multiple constructs and in multiple expression systems. Also, when optimizing the protein solubility and stability, they need to test many different conditions. In traditional small scale high throughput purification technologies researchers often face challenges with reproducibility, purity and recovery. Here we show examples where PhyTip® columns overcome some of these challenges.

Figure 4 displays some data from a study performed at Roche. This group consistently clones and expresses varying constructs as part of an X-ray crystallography screen and must test different expression systems to identify forms with the highest expression levels, stability, and potency. In this study they were looking at the incorporation of PhyTip® columns in their process. The SDS-PAGE gels demonstrate a level of purity and consistency of the purified constructs across their entire screen.

A. 16 Constructs expressed in E. coli cultures.



B. 9 Constructs expressed in Baculoviral cultures.



C. 3 Constructs expressed in Mammalian cultures.




Figure 4. SDS-PAGE gels on protein purity from different constructs expressed in different expression systems.

A.

† pH HEPES/
NaOH
 ‡ mM HEPES/
NaOH
 • %v/v Glycerol
 † mM NaCl
 † pH Tris-HCl/
NaOH
 † mM Tris-HCl/
NaOH
 † mM L-Arginine


A1	A2	A3	A4	A5	A6	A7	A8	A9	A10	A11	A12
7.00 † 20.0 ‡ 5.0 • 150.0 †	7.00 † 20.0 ‡ 5.0 • 325.0 †	7.00 † 20.0 ‡ 5.0 • 500.0 †	7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †
B1	B2	B3	B4	B5	B6	B7	B8	B9	B10	B11	B12
7.00 † 20.0 ‡ 5.0 • 150.0 †	7.00 † 20.0 ‡ 5.0 • 325.0 †	7.00 † 20.0 ‡ 5.0 • 500.0 †	7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †
C1	C2	C3	C4	C5	C6	C7	C8	C9	C10	C11	C12
7.00 † 20.0 ‡ 5.0 • 150.0 †	7.00 † 20.0 ‡ 5.0 • 325.0 †	7.00 † 20.0 ‡ 5.0 • 500.0 †	7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †
D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	D11	D12
7.00 † 20.0 ‡ 5.0 • 150.0 †	7.00 † 20.0 ‡ 5.0 • 325.0 †	7.00 † 20.0 ‡ 5.0 • 500.0 †	7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †
E1	E2	E3	E4	E5	E6	E7	E8	E9	E10	E11	E12
7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †	8.50 † 20.0 ‡ 5.0 • 150.0 †	8.50 † 20.0 ‡ 5.0 • 325.0 †	8.50 † 20.0 ‡ 5.0 • 500.0 †
F1	F2	F3	F4	F5	F6	F7	F8	F9	F10	F11	F12
7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †	8.50 † 20.0 ‡ 5.0 • 150.0 †	8.50 † 20.0 ‡ 5.0 • 325.0 †	8.50 † 20.0 ‡ 5.0 • 500.0 †
G1	G2	G3	G4	G5	G6	G7	G8	G9	G10	G11	G12
7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †	8.50 † 20.0 ‡ 5.0 • 150.0 †	8.50 † 20.0 ‡ 5.0 • 325.0 †	8.50 † 20.0 ‡ 5.0 • 500.0 †
H1	H2	H3	H4	H5	H6	H7	H8	H9	H10	H11	H12
7.50 † 20.0 ‡ 5.0 • 150.0 †	7.50 † 20.0 ‡ 5.0 • 325.0 †	7.50 † 20.0 ‡ 5.0 • 500.0 †	7.80 † 20.0 ‡ 5.0 • 150.0 †	7.80 † 20.0 ‡ 5.0 • 325.0 †	7.80 † 20.0 ‡ 5.0 • 500.0 †	8.00 † 20.0 ‡ 5.0 • 150.0 †	8.00 † 20.0 ‡ 5.0 • 325.0 †	8.00 † 20.0 ‡ 5.0 • 500.0 †	8.50 † 20.0 ‡ 5.0 • 150.0 †	8.50 † 20.0 ‡ 5.0 • 325.0 †	8.50 † 20.0 ‡ 5.0 • 500.0 †

Figure 5a. Buffer composition of samples for testing protein stability and binding efficiency.

Tris- HCl pH8.5 buffer containing different concentrations of NaCl and detergents




Hepes, pH7.5 buffers containing different concentrations of NaCl and detergents



a = 150mM NaCl; b = 300mM NaCl; c = 500mM NaCl


1 = 0.1% Igepal; 2 = 0.1% n-Octyl-β-D-glucopyranoside; 3 = 0.5% Igepal; 4 = 0.2% n-Dodecyl-β-D-maltoside

Figure 5b. Impact of NaCl concentrations and detergents.**Figure 5c.** Impact of arginine for binding efficiency.

As a conclusion this application example showed how PhyTip® columns can be used to enable high throughput purification to evaluate binding and stability characteristics with high levels of consistency and purity.

Immunoprecipitation (IP) is a widely used technique to isolate specific types of proteins or a part of a protein containing a specific biomarker from cell or tissue lysates. It is essential for studying post-translational modifications (PTMs) as well as studying drug-target interactions. A specific antibody is immobilized on a bead matrix and used to bind its antigen. Traditional immunoprecipitation is performed in a microcentrifuge tube and each step is carried out in a batch-binding manner with the use of a centrifuge.

Translating the manual IP process to the PhyTip® column is highlighted here by the work lead by Lilian Phu at Genentech.






Figure 7. The immunoprecipitation method was used to isolate ubiquitin remnant motif (KGG) peptides from a complex mixture. These specific peptides were then analyzed by LC-MS/MS for KGG motif identity and for the level of non-specific contaminating peptides. When comparing the automated PhyTip® column process with the manual process, the automated procedure was more comparable. The automated process displayed more consistency and a lower background of non-specific peptides.

Conclusion

PhyTip® columns are a versatile platform that can be adapted to many different biomolecule purification applications. They can be used on all major automated liquid handling systems, enabling efficient set up of purification runs. Compared to alternative technologies PhyTip® columns offer easy automation, high sample recoveries and optimized concentrations through small elution volumes. This allows for high throughput solutions and purified biomolecule samples optimally prepared for following characterization steps.

Additional Application Data on PhyTip® Columns

Application Note: “MiniPrep PhyTip® Columns for Plasmid DNA Purification”

Application Note: “ProPlus PhyTip® Affinity Columns... The Antibody Purification System that is the Best of Protein A and Protein G”

Application Note: “Automated High Throughput Desalting and Buffer Exchange PhyTip® Columns”

Application Note: True parallel purification using PhyTip® column and Tecan Freedom EVO with Multichannel arm 96

Application Note: “High Performance Immunoprecipitation (HPIP)” Direct method

Application Note: “High Performance Immunoprecipitation (HPIP)” Indirect method

Application Note: “Optimization Strategies for High Performance Purification and Analysis of Recombinant Proteins with Micro Volume PhyTip® Columns and Caliper Life Sciences Automation”

Application Note: “Perform Resin Screening for Most Efficient Capture of Target Proteins”

EUROPE

Main Office: +46 18 565900
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 800 446 4752
Order Fax: +1 704 654 4917
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: +82 31 706 8500
Fax: +82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

INDIA

Tel: +91 11 45653772
india@biotage.com

Distributors in other regions
are listed on www.biotage.com

Literature Number: PPS696

© 2022 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.