# Biotage® PhyPrep

Maxi, Mega, Giga





# **Automated Plasmid Purification**

Had enough of manual, time-consuming plasmid purifications? So have we. So we invented the Biotage® PhyPrep, which automates the purification of plasmid DNA from cell pellets at the Maxi, Mega and Giga scales. The result is highly pure, endotoxin-free\* plasmid DNA – the optimal starting point for downstream applications, such as transient transfection. Biotage® PhyPrep comes with everything needed to generate highly pure plasmids, from a simple sample preparation step to a fully automated purification procedure, freeing up your time to get on with the really important work – research.



- ✓ Transfection ready plasmid DNA
- ✓ Hands-free walkaway automation
- ✓ Intuitive, hassle-free software and workflow
- ✓ Consistent quality and results

Plasmid yield

1 mg (maxiprep) 5 mg (megaprep) 10 mg (gigaprep)

||||

Samples per run 1-4

••

Endotoxin levels <0.1 EU/μg\*



PhyPrep takes care of the plasmids, you can focus on your science.


\*With rich media and high copy plasmids

# **Endotoxin-free Transfection Grade Plasmid DNA**

We know that when it comes to plasmid quality, there's no room for compromise. Using our proprietary dual flow chromatography technology, Biotage® PhyPrep provides endotoxin-free\*, supercoiled plasmid DNA. This highly pure plasmid DNA is suitable for many downstream applications, such as transient transfection, in vivo studies, cloning, genome editing and other molecular biology assays.

\*With rich media and high copy plasmids

## **Maxiprep Results**



Lane	Sample
1	DNA Ladder
2	Normalized Qiagen
3	Normalized Maxiprep #1
4	Normalized Maxiprep #2
5	Normalized Maxiprep #3
6	Normalized Maxiprep #4


Sample ID	Yield (mg)	Concentration (ng/uL)	A260/ A280	A260/ A230	Volume (mL)	Endotoxin Level (EU/µg)
Maxiprep #1	1.46	301	1.90	2.31	4.85	< 0.02
Maxiprep #2	1.66	334	1.91	2.31	4.95	< 0.02
Maxiprep #3	1.55	322	1.90	2.31	4.80	< 0.02
Maxiprep #4	1.60	287	1.92	2.31	5.57	< 0.02

**Table 1.** UV absorbance data of plasmids purified with the PhyPrep MaxiPrep Kit Actual yield is dependent on the plasmid copy number, culture growth conditions, strain of *E.coli* utilized, and cell pellet weight processed.

Agarose gel electrophoresis analysis of plasmid DNA purified by the PhyPrep MaxiPrep Kit.

The gel was run on a 1% Agarose gel at 70 volts for 2.5 hours, then stained for 10 minutes and destained for 20 minutes. All samples were normalized.

# Megaprep Results



Agarose gel electrophoresis analysis of plasmid DNA
purified by the PhyPrep MegaPrep Kit.


The gel was run on a 1% Agarose gel at 70 volts for 2.5 hours, then stained for 10 minutes and destained for 20 minutes. All samples were normalized.

Lane	Sample
1	DNA Ladder
2	Normalized Qiagen
3	Normalized Megaprep #1
4	Normalized Megaprep #2
5	Normalized Megaprep #3
6	Normalized Megaprep #4

Sample ID	Yield (mg)	Concentration (ng/uL)	A260/ A280	A260/ A230	Volume (mL)	Endotoxin Level (EU/µg)
Megaprep #1	5.91	326	1.92	2.42	18.1	< 0.02
Megaprep #2	5.9	330	1.92	2.43	17.9	< 0.02
Megaprep #3	5.0	270	1.91	2.39	18.6	< 0.02
Megaprep #4	6.0	301	1.91	2.35	20	< 0.02

**Table 2.** UV absorbance data of plasmids purified with the PhyPrep MegaPrep Kit Actual yield is dependent on the plasmid copy number, culture growth conditions, strain of *E.coli* utilized, and cell pellet weight processed.

# Gigaprep Results



Agarose gel electrophoresis analysis of plasmid DNA purified by the PhyPrep GigaPrep Kit.


The gel was run on a 1% Agarose gel at 70 volts for 2.5 hours, then stained for 10 minutes and destained for 20 minutes. All samples were normalized.


Lane	Sample
1	DNA Ladder
2	Normalized Qiagen
3	Normalized Gigaprep #1
4	Normalized Gigaprep #2
5	Normalized Gigaprep #3
6	Normalized Gigaprep #4


Sample ID	Yield (mg)	Concentration (ng/uL)	A260/ A280	A260/ A230	Volume (mL)	Endotoxin Level (EU/µg)
Gigaprep #1	16.55	537	1.98	2.38	30.82	0.009
Gigaprep #2	11.94	379	1.94	2.34	31.53	0.070
Gigaprep #3	10.41	363	1.92	2.35	28.71	0.014
Gigaprep #4	13.75	530	2.01	2.41	25.96	0.030

**Table 3.** UV absorbance data of plasmids purified with the PhyPrep GigaPrep Kit Actual yield is dependent on the plasmid copy number, culture growth conditions, strain of *E.coli* utilized, and cell pellet weight processed.

## **Transient Transfection Results**








**Results from transfection experiment performed by Mirus Bio.** Successful transfection is determined via the percent of cells positive for GFP. In all cases, plasmids successfully transfected over 95% of the mammalian cells. pEFGP is an external plasmid used as a standard control for the transfection reagent.















# Simple and intuitive


# **Automation Made Easy**

Biotage® PhyPrep is designed to make automation simple and easy to introduce to your laboratory. Cell pellet preparation makes use of pre-aliquoted buffers so that you can get your samples prepared within minutes in a reproducible manner without the need to dispense and measure solutions. The software is incredibly easy to use, requiring only a few options to select on the touch screen to set your purifications in motion. And once the PhyPrep is dealing with your sample, you can leave it, confident that when you return you'll have endotoxin-free\* purified plasmids.









\*With rich media and high copy plasmids


# **Inspired Chromatography**

## The Heart of the Process

Automation provides a simple and reliable way to purify plasmids, taking a laborious and manual process away from highly-skilled scientists and freeing up time for more valuable work. This is achieved through the clever application of chromatographic purification, the heart of the process.

## Cutting-edge Chromatography

Automated plasmid purification is made possible through our proprietary Biotage® Växel columns, packed with anion exchange resin. The sample is processed by dual flow chromatography — in this unique cutting-edge process, the sample flows back and forth through the Växel column. Dual flow chromatography is a purification process that we have employed for small-scale automated protein purification since 2004.



# **Technical Specifications**

### **Operational Performance Specifications**

Purification scale Maxipre	o, Megaprep and Gigaprep
Plasmid Vield	5/10 mg from 5/8/16 g cell wet weight using rich media and a high-copy at Maxi/Mega/Giga scale
Throughput 1–4 sam	ples per run

Ш	Throughput	1-4 samples per ru	n				
		Prep Scale	Elution Volume	Number of Sample Preps	Run Time (Hours/Minutes)		
				1	42 m		
		Massinson	5 mL	2	42 m		
		Maxiprep		3	1 h 10 m		
				4	1 h 10 m		
<b>₹</b> \	Processing time			1	1 h 16 m		
	Frocessing time	Megaprep	18 mL	2	1 h 16 m		
		медаргер	10 IIIL	3	2 h 22 m		
				4	2 h 22 m		
				1	2 h 31 m		
		Gigaprep	28 ml	2	2 h 34 m		
		Gigapi ep	28 mL	3	4 h 57 m		
				4	4 h 58 m		
I	Elution volume	5/18/28 mL for Maxiprep, Megaprep, and Gigaprep					
•	Endotoxin levels of purified plasmid DNA	<0.1 EU/µg (with rich media and high copy plasmids)					
	Purity of purified plasmid DNA	A260/A280 > 1.8, A260/A230 > 2.0					
enei	ral System Specifications						
	Weight	23 kg (50 lbs)					
<b>.</b>	Dimensions	Footprint: (W x D) 530 mm x 710 mm (20.8" x 27.9") Height: 830 mm (32.6")					
Ł)	Interfaces	Touch screen 12"					
4	Max Power Consumption	580 VA					
Ç	Certifications	CE, CB, TUV, NRTL					

# **Ordering Information**

	5	
Part number	Description	Qty
Instrument		
416700	Biotage PhyPrep (Requires 356330SP)	1
356330SP	Vacuum pump ME1C	1
Consumables		
DPM-16-05-72-KIT	Biotage Plasmid MaxiPrep kit Endotoxin- Free (16 Samples)	1
DPM-16-10-72-KIT	Biotage Plasmid MegaPrep kit Endotoxin- Free (16 Samples)	1
DPM-08-20-72-KIT	Biotage Plasmid GigaPrep kit Endotoxin- Free (8 samples)	1
417458SP	50 x PhyPrep Tubes for Biotage® PhyPrep	1
DPM-00-10-00	$100 \times 20 \text{ mL transfer tips for Biotage}^{\circ}$ PhyPrep	1
DPM-00-25-00	25 x 20 mL transfer tips for Biotage® PhyPrep	1



# Your Complete Partner for Automated Protein and Plasmid Purification

Biotage is a worldwide supplier of instruments and accessories designed to facilitate the work of scientists in life sciences. With our deep knowledge of the industry, academic contacts and in-house R&D teams, we can deliver the best solutions to your challenges. We take great pride in our flexibility and ability to meet our customer's individual needs. With strong foundations in analytical, organic, process, and biomolecule chemistry, we can offer the widest range of solutions available on the market.

#### EUROPE

Main Office: +46 18 565900
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

#### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 800 446 4752
Order Fax: +1 704 654 4917
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
us-1-pointsupport@biotage.com

#### JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp\_order@biotage.com
jp-1-pointsupport@biotage.com

# **CHINA**Tel: +86 21 68162810

Fax: +86 21 68162829 cn\_order@biotage.com cn-1-pointsupport@biotage.com

#### KOREA

Tel: +82 31 706 8500 Fax: +82 31 706 8510 korea\_info@biotage.com kr-1-pointsupport@biotage.com

#### INDIA

Tel: +91 11 45653772 india@biotage.com

Distributors in other regions are listed on www.biotage.com

