

Biotage® Solvent Recovery SVOC



The Biotage® Solvent Recovery SVOC is designed to condense and collect solvent vapors which are generated when using the DryVap® concentration system.

Biotage® Solvent Recovery SVOC can help recover solvents evaporated during the concentration step prior to injection into the GC or HPLC. The Solvent Recovery SVOC comes coupled with the DryVap® in-line drying and evaporation system for a complete dry-evaporate-trap solvent system, consistent with the workflow and environmentally friendly.

The system can recover 95% dichloromethane(DCM) evaporated from six 180 mL samples simultaneously, using vacuum and an auxilliary chiller.

Advantages

- » Typical solvent volumes of 180–200 mL of DCM can be evaporated in 60 minute run times
- » The end point can be automatically detected and the process ended when the endpoint is reached
- » The sample can be evaporated directly into GC autosampler vials, eliminating transfer from the evaporation vessel
- » Excellent solvent recovery for disposal or re-use

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: + 46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: + 82 31 706 8500
Fax: + 82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

INDIA

Tel: +91 22 4005 3712
india@biotage.com

Distributors in other regions are listed on www.biotage.com

Literature Number: PPS532

© 2019 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.