

Sample Preparation Guide

Required Kit



RNase A



Resuspension Buffer



Lysis Buffer




Precipitation Buffer




Thawed Cell Pellet

1 Resuspend Cells




2 Lyse Cells



3 Precipitate



4 Filter



Literature Number: 417480-A

© 2021 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.