

# Sole Source Document

## Biotage® Extrahera™ LV-200

28 April 2021

Biotage is the sole manufacturer and supplier of the Extrahera™ LV-200 automation system for solid phase extraction, supported liquid extraction, phospholipid depletion, protein precipitation, and dual flow chromatography with the following specifications.

Patent is pending for the design of the collection plate lift which is used to eliminate/minimise potential cross-contamination.

### Unique features of the Extrahera™ LV-200 include:

- » Two tier design with lower carousel to minimise footprint.
- » Integral fume hood capability to allow bench top installation (with appropriate exhaust)
- » Self-refilling solvent reservoirs
- » Dual flow positive pressure processing for columns and plates
- » Dual flow chromatography for in-tip sample prep
- » Re-use of solvent pipette tips where appropriate
- » Ability to use multiple tip configuration on processing shelf for SLE and SPE



#### EUROPE

Main Office: +46 18 565900  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: +46 18 56 57 11  
 eu-1-pointsupport@biotage.com

#### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 800 446 4752  
 Order Fax: +1 704 654 4917  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 us-1-pointsupport@biotage.com

#### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

#### CHINA

Tel: +86 21 68162810  
 Fax: +86 21 68162829  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

#### KOREA

Tel: +82 31 706 8500  
 Fax: +82 31 706 8510  
 korea\_info@biotage.com  
 kr-1-pointsupport@biotage.com

#### INDIA

Tel: +91 11 45653772  
 india@biotage.com

Distributors in other regions are listed on [www.biotage.com](http://www.biotage.com)

#### Literature Number: PPS648

© 2021 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.